

| Shade VITA® Classical sorting by color groups → | | A0,5 | | | A1 | | | A2 | | | A3 | | | A3,5 | | | A4 | | | B1 | | | B2 | | | B3 | | | B4 | | | C1 | | | C2 | | | C3 | | | C4 | | | D2 | | | D3 | | | D4 | | | | | |
|---|----|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|
| Thoot Area → | | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal |
| ASTERIA shade ↓ | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A1B | VG | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A2B | | VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A3B | | | | VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A3,5B | | | | | | VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A4B | | | | | | VG | | VG VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| BL | VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| B3B | | | | | | | | | | | | G VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| NE | | G | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | | | | | | | | |
| OcE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| WE | | VG | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | | | | | | | | | | |
| TE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| YE | | | | | | G | | G | | G | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

| Shade VITA® Classical sorting by brightness (value) → | | A0,5 | | | B1 | | | A1 | | | B2 | | | D2 | | | A2 | | | C1 | | | C2 | | | D4 | | | A3 | | | D3 | | | B3 | | | A3,5 | | | B4 | | | C3 | | | A4 | | | C4 | | | | | |
|---|----|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|-------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|----------|------|---------|
| Thoot Area → | | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal | Cervical | Body | Incisal |
| ASTERIA shade ↓ | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A1B | VG | G VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A2B | | VG | | VG | | G VG | | G VG | | G VG | | G VG | | G VG | | G VG | | G VG | | G VG | | G VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| A3B | | | | | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | G G | | G VG | | G VG | | G G | | | | | | | | | | | | | | | | | | | | | | | | | |
| A3,5B | | | | | | | | | | | | | | | | | | | | | | | | G | | VG | | VG | | VG | | G VG | | G G | | G | | | | | | | | G VG | | | | | | | | | | | |
| A4B | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | VG | | VG | | | | | | VG VG | | VG | | | | | | | | | | | | | | | |
| BL | VG | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| B3B | | | | | | | | | | | | | | | | | | | | | | | | VG VG | | | | | | G VG | | G VG | | VG VG | | | | | | | | | | | | | | | | | | | | | |
| NE | | G | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | VG | | | | | | | | | |
| OcE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| WE | | VG | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | | | | | | | | | | |
| TE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| YE | | | | | | | | | | | | | | | | | | | | | | | | | | G | | G | | G | | G | | G | | G | | G | | G | | G | | G | | | | | | | | | | | |